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
12 REFERENZEN / LITERATUR

1. Jopling H, Yates A, Burgoyne N, Hayden K, Chaloner C, Tetlow L. Paediatric Anti-Müllerian Hormone measurement: Male and female reference intervals established using the automated Beckman Coulter Access AMH assay. *Endocrinol Diabetes Metab.* 2018;1(4):e00021.
2. Hampl R, Šnajderová, M, Mardešić T. Antimüllerian Hormone (AMH) Not Only a Marker for Prediction of Ovarian Reserve. *Physiol. Res.* 60: 217-223, 2011.
3. Matuszczak E, Hermanowicz A, Komarowska M, Debek W. Serum AMH in Physiology and Pathology of Male Gonads. *Int J Endocrinol.* 2013:128907.
4. Victoria M, Labrosse J, Krief F, Cédric-Durnerin I, Comtet M, Grynberg M. Anti Müllerian Hormone: More than a biomarker of female reproductive function. *J Gynecol Obstet Hum Reprod.* 2019 Jan;48(1):19-24.
5. Broer S, Broekmans F, Laven J, Fauser B. Anti-Müllerian hormone: ovarian reserve testing and its potential clinical implications. *Hum Reprod Update.* 2014 Sep-Oct;20(5):688-701.
6. Broer S, Eijkemans M, Scheffer G, van Rooij I, de Vet A, Themmen A, Laven J, de Jong F, te Velde E, Fauser B, Broekmans F. Anti-Müllerian Hormone Predicts Menopause: A Long-Term Follow-Up Study in Normoovulatory Women. *J Clin Endocrinol Metab.* 2011 Aug;96(8):2532-9.
7. Visser J, de Jong F, Laven J, Themmen A. Anti-Müllerian hormone: a new marker for ovarian function. *Reproduction.* 2006 Jan;131(1):1-9.
8. Broer S, van Disseldorp J, Broeze K, Dolleman M, Opmeer B, Bossuyt P, Eijkemans M, Mol B, Broekmans F; IMPORT study group. Added value of ovarian reserve testing on patient characteristics in the prediction of ovarian response and ongoing pregnancy: an individual patient data approach. *Hum Reprod Update.* 2013 Jan-Feb;19(1):26-36.
9. Ebner T, Sommergruber M, Moser M, Shebl O, Schreier-Lechner E, Tews G. Basal level of anti-Müllerian hormone is associated with oocyte quality in stimulated cycles. *Hum Reprod.* 2006 Aug;21(8):2022-6.
10. Sun H, Mao H, Cai J, Zhao Y. Research progress on anti-müllerian hormone clinical applications and immunoassay development. *Frontiers in Laboratory Medicine.* 2018 2(1):14-18.
11. Anckaert E, Öktem M, Thies A, Cohen-Bacrie M, Daan NM, Schiettecatte J, Müller C, Topcu D, Gröning A, Ternaux F, Engel C, Engelmann S, Milczynski C. Multicenter analytical performance evaluation of a fully automated anti-Müllerian hormone assay and reference interval determination. *Clinical biochemistry* 2016, 49(3), 260–267.
12. Anderson R A, Anckaert E, Bosch E, Dewailly D, Dunlop CE, Fehr D, Nardo L, Smitz J, Tremellen K, Denk B, Geistanger A, Hund M. Prospective study into the value of the automated Elecsys antimüllerian hormone assay for the assessment of the ovarian growing follicle pool. *Fertility and sterility* 2015, 103(4), 1074–1080.e4.
13. Practice Committee of the American Society for Reproductive Medicine. Testing and interpreting measures of ovarian reserve: a committee opinion. *Fertility and sterility* 2012, 98(6), 1407–1415.

Please use only the valid version of the Instructions for Use provided with the kit

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