



# CERTIFICATION

**AOAC<sup>®</sup> Performance Tested<sup>SM</sup>**

Certificate No.

**021402**

The AOAC Research Institute hereby certifies that the performance of the test kit known as:

**HistaSure<sup>™</sup> ELISA<sup>Fast Track</sup>**

manufactured by

**LDN Labor Diagnostika Nord GmbH & Co  
KG Am Eichenhain 1  
48531 Nordhorn  
Germany**

This method has been evaluated in the AOAC<sup>®</sup> *Performance Tested Methods<sup>SM</sup>* Program, and found to perform as stated by the manufacturer contingent to the comments contained in the manuscript. This certificate means that an AOAC<sup>®</sup> Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC *Performance Tested<sup>SM</sup>* certification mark along with the statement - "THIS METHOD'S PERFORMANCE WAS REVIEWED BY AOAC RESEARCH INSTITUTE AND WAS FOUND TO PERFORM TO THE MANUFACTURER'S SPECIFICATIONS" - on the above mentioned method for a period of one calendar year from the date of this certificate (November 05, 2018 – December 31, 2019). Renewal may be granted at the end of one year under the rules stated in the licensing agreement.

*Scott Coates*

\_\_\_\_\_  
Scott Coates, Senior Director  
Signature for AOAC Research Institute

\_\_\_\_\_  
November 05, 2018

Date

**METHOD AUTHORS**

Georg Manz and Esther Bootlink

**SUBMITTING COMPANY**LDN Labor Diagnostika Nord GmbH & Co. KG  
Am Eichenhain 1  
48531 Nordhorn  
Germany**KIT NAME(S)**HistaSure™ ELISA<sup>Fast Track</sup>**CATALOG NUMBERS**

FC E-3600, FC E-3900

**INDEPENDENT LABORATORY**Eurofins Central Analytical Laboratories  
2315 N. Causeway Blvd.  
Metairie, LA 70001  
USA**AOAC EXPERTS AND PEER REVIEWERS**Christina A. Mireles DeWitt<sup>1</sup>, Sevim Köse<sup>2</sup>, James Hungerford<sup>3</sup>  
<sup>1</sup> Oregon State University, Astoria, OR USA  
<sup>2</sup> Karadeniz Technical University, Trabzon, Turkey  
<sup>3</sup> U.S. Food and Drug Administration, Bothel, WA, USA**APPLICABILITY OF METHOD**

Target analyte– histamine

Matrices – (10 g samples) - fresh/frozen yellowfin tuna, canned tuna-chunk light in water, fresh/frozen mahi mahi, canned sardines in oil, and fish meal

**Performance claims**

Selectivity: Cross reactivity with components of similar behavior &lt; 2%

Accuracy: %RSD &lt; 15% (above LOQ)

Precision: recovery within 80 – 120 %

Detection limit: LOQ &lt; 3 ppm (Corresponding to Control 2 of the ELISA)

Stability: 2 year shelf life at 2 – 8°C. Acceptable OD range of control 0: 1.250 – 2.750.

Method comparison LDN HistaSure™ ELISA<sup>Fast Track</sup> vs AOAC 977.13.Correlation Coefficient r>0.9. Ratio AOAC 977.13 / HistaSure™ ELISA<sup>Fast Track</sup> 0.8 – 1.2**REFERENCE METHOD**W. Horwitz, G.W. Latimer, 2005. AOAC Official Method 977.13. Official Methods of analysis of AOAC International 18<sup>th</sup> Edition 35.1.32 (13)**ORIGINAL CERTIFICATION DATE**

February 03, 2014

**CERTIFICATION RENEWAL RECORD**

Renewed annually through 2019

**METHOD MODIFICATION RECORD**

1. November 2017 Level 1

**SUMMARY OF MODIFICATION**

1. Editorial changes

Under this AOAC® *Performance Tested*<sup>SM</sup> License Number, 021402 this method is distributed by:  
NONEUnder this AOAC® *Performance Tested*<sup>SM</sup> License Number, 021402 this method is distributed as:  
NONE**PRINCIPLE OF THE METHOD (1)**

Histamine is extracted from the fish using a simple water extraction. After filtration or centrifugation the clear supernatant is used for the subsequent step, the derivatization of histamine.

By use of the acylation reagent, histamine from the calibrators and histamine from the food extracts is quantitatively derivatized into N-acylhistamine. The competitive Histamine ELISA kit uses the microtiter plate format. Histamine is bound to the solid phase of the microtiter plate. Acylated histamine and solid phase bound histamine compete for a fixed number of antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antiserum-peroxidase complexes are removed by washing. The substrate TMB/peroxidase reaction is monitored at 450 nm. The amount of antibody bound to the solid phase histamine is inversely proportional to the histamine concentration of the sample.

**DISCUSSION OF THE VALIDATION STUDY (1)**In the AOAC PTM validation study the performance of the LDN HistaSure™ ELISA<sup>Fast Track</sup> was at least equivalent to the AOAC 977.13 reference method. The LDN HistaSure™ ELISA<sup>Fast Track</sup> showed good recovery rates for all tested matrices (fresh/frozen tuna, canned tuna, fresh/frozen mahi mahi, canned sardines and fish meal) and repeatability precision was excellent at the 50 ppm defect action level. No solvents are required for extraction. The method is robust, fast, efficient and suitable for high throughput screening. Based on the results of this study it could be recommended that the LDN HistaSure™ ELISA<sup>Fast Track</sup> receives PTM certification.

**Table 9: Method comparison results overview AOAC 977.13 vs. LDN HistaSure™ ELISA <sup>Fast Track</sup> (1)**

	Fortification level (ppm)	AOAC 977.13				LDN HistaSure™ ELISA <sup>Fast Track</sup>			
		Mean concentration (n=5) (ppm)	SD	%RSD	Average % recovery (n=5)	Mean concentration (n=7) (ppm)	SD	%RSD	Average % recovery (n=7)
Fresh/frozen Tuna	0	0.51	0.63	123.5	/	0.27 (< LOQ)	0.11	38.7	/
	5	4.92	0.81	16.4	88.1	5.25	0.40	7.52	99.6
	20	17.6	2.81	16.0	85.4	19.2	2.08	10.8	94.6
	50	43.9	7.37	16.8	86.8	42.9	1.77	4.13	85.2
	100	87.9	15.1	17.2	87.4	92.6	6.99	7.55	92.4
	250	261.8	27.2	10.4	104.5	218.6	17.3	7.90	87.3
	<b>Mean recovery</b>				<b>90.4</b>				<b>91.8</b>
Canned tuna	0	0.22	0.22	99.7	/	0.50 (< LOQ)	0.17	33.7	/
	5	4.84	0.16	3.21	92.4	5.49	0.76	13.9	99.8
	20	18.3	2.35	12.8	90.6	20.2	2.01	9.92	98.6
	50	52.3	2.89	5.52	104.2	47.8	4.10	8.57	94.7
	100	99.5	7.12	7.15	99.3	99.8	8.10	8.12	98.7
	250	238.9	33.8	14.2	95.5	267.3	8.46	3.17	106.5
	<b>Mean recovery</b>				<b>96.4</b>				<b>99.8</b>
Frozen mahi mahi	0	1.52	0.96	63.3	/	1.22 (< LOQ)	0.21	17.3	/
	5	5.44	0.77	14.3	78.4	6.38	0.61	9.56	103.1
	20	17.9	1.69	9.45	82.0	18.9	1.28	6.78	88.2
	50	45.8	7.68	16.8	88.5	43.5	2.29	5.26	84.6
	100	93.9	6.30	6.71	92.4	82.7	2.46	2.97	81.5
	250	242.6	19.9	8.22	96.4	199.0	12.6	6.35	79.1
	<b>Mean recovery</b>				<b>87.5</b>				<b>87.3</b>
Canned sardines	0	0.21	0.32	155.1	/	0.67 (< LOQ)	0.19	28.9	/
	5	4.24	0.28	6.57	80.6	5.27	0.28	5.36	92.1
	20	20.7	1.10	5.32	102.3	18.2	1.24	6.78	87.8
	50	41.6	4.89	11.8	82.9	41.5	0.75	1.80	81.6
	100	86.1	10.8	12.5	85.9	77.5	3.94	5.08	76.8
	250	230.3	32.2	14.0	92.0	249.9	23.0	9.22	99.7
	<b>Mean recovery</b>				<b>89.0</b>				<b>87.6</b>
Fish meal	0	6.31	0.73	11.6	/	9.4	0.54	5.71	/
	5	5.30	0.30	5.73	-20.1	14.0	1.11	7.92	92.3
	20	15.5	1.42	9.12	46.2	25.5	1.02	3.98	80.5
	50	38.5	4.11	10.7	64.5	48.9	1.02	2.09	79.0
	100	95.1	2.78	2.92	88.8	93.8	3.73	3.97	84.4
	250	212.9	19.8	9.29	82.6	244.2	15.8	6.47	93.9
	<b>Mean recovery</b>				<b>52.4</b>				<b>86.0</b>
Naturally contaminated	Sample 1	61	3.06	5.04	/	58	5.2	8.95	/
	Sample 2	203	13.7	6.74	/	133	10.7	7.99	/
	Sample 3	7	1.00	14.3	/	8.6	0.92	10.7	/
	Sample 4	54	3.21	5.99	/	48	2.8	5.83	/
	Sample 5	57	1.53	2.66	/	69	5.3	7.75	/

#### REFERENCES CITED

1. Manz, Georg, and Boolsink, Esther., Evaluation of the HistaSure™ ELISA <sup>Fast Track</sup> for Determination of Histamine in Fish Samples, AOAC® *Performance Tested*<sup>SM</sup> certification number 021402.
2. AOAC Research Institute Validation Outline for HistaSure™ ELISA <sup>Fast Track</sup>, Approved –February 2014.
3. P.L. Rogers, W.F. Staruszkiewicz, 2000. Histamine Test Kit Comparison. *Journal of Aquatic Food Product Technology* 9 (2), 5-17
4. J.M. Hungerford, 2010. Scombroid poisoning: a review. *Toxicon* 15, 231-243
5. S. Köse, N. Kaklikkaya, S. Koral, B. Tufan, K.C. Buruk, F. Aydin, 2011. Commercial test kits and the determination of histamine in traditional (ethnic) fish products- evaluation against EU accepted HPLC method. *Food Chemistry*, 125 (4) 1490-1497
6. A. Lupo, M. Mozola, 2011. Validation study of a rapid ELISA for the detection of histamine in tuna. *Journal of AOAC International* 94 (3), 886-899
7. D.G. Altman. J.M. Bland, 1983. Measurement in medicine: the analysis of method comparison studies. *The Statistician* 32, 307-317
8. R. Taylor. 1990. Interpretation of the correlation coefficient: A basic review. *Journal of diagnostic medical sonography* 6(1), 35-39
9. National Association of Testing Authorities, Australia, 2012. Technical Note 17: Guidelines for the validation and verification of quantitative and qualitative test methods.
10. US department of health and human services food and drug administration, 2001. Guidance for Industry, bioanalytical method validation.
11. M. Thompson, S.L.R. Ellison, R. Wood, 2002. Harmonized guidelines for single laboratory validation of methods analysis. *Pure and applied chemistry* 74, 835-855
12. NCCLS, Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. EP 17-A Vol. 24 No. 10
13. W. Horwitz, G.W. Latimer, 2005. AOAC Official Method 977.13. Official Methods of analysis of AOAC International 18<sup>th</sup> Edition 35.1.32